

	Dialyzer	Surface Area (m2)	Priming Volume (ml)	Clearance Qb300, Qd500				Clearance Qb400, Qd500				Clearance Qb400, Qd800				KoA Urea
				Urea	Creatinine	Phosphate	Vit B12	Urea	Creatinine	Phosphate	Vit B12	Urea	Creatinine	Phosphate	Vit B12	
Dimesol	B-16H	1.6	88	267	250	235	175	320	295	265	190	–	–	–	–	1064.9
Fresenius	Optiflux F160NR	1.5	83	266	238	230	152	308	264	275	157	344	295	296	173	1064
Fresenius	Optiflux F180NR	1.8	99	274	251	238	168	323	286	287	178	358	316	314	203	1239
Gambro	Polyflux 140H	1.4	94	262	232	220	149	309	266	250	163	–	–	–	–	1040
Gambro	Polyflux 170H	1.7	115	270	243	232	162	321	281	266	178	–	–	–	–	1190
Baxter	Exeltra 150	1.5	95	262	242	227	152	305	274	255	163	–	–	–	–	993
Baxter	Exeltra 170	1.7	105	268	252	232	160	310	286	261	172	–	–	–	–	1103

	Dialyzer	Surface Area (m2)	Priming Volume (ml)	Urea	Creatinine	Phosphate	Vit B12	Urea	Creatinine	Phosphate	Vit B12	Urea	Creatinine	Phosphate	Vit B12	KoA Urea
Baxter	Exeltra 170	1.7	105	268	252	232	160	310	286	261	172	–	–	–	–	1103
Fresenius	Optiflux F180NR	1.8	99	274	251	238	168	323	286	287	178	358	316	314	203	1239
Baxter	Exeltra 190	1.9	115	273	251	242	168	323	289	276	183	–	–	–	–	1214
Fresenius	Optiflux F200NR	2.0	112	277	253	250	173	330	289	290	189	364	325	319	206	1317
Gambro	Polyflux 210H	2.1	125	281	259	249	183	339	303	289	203	–	–	–	–	1500

	Dialyzer	Surface Area (m2)	Priming Volume (ml)	Urea	Creatinine	Phosphate	Vit B12	Urea	Creatinine	Phosphate	Vit B12	Urea	Creatinine	Phosphate	Vit B12	KoA Urea
Fresenius	Optiflux F200NR	2.0	112	277	253	250	173	330	289	290	189	364	325	319	206	1317
Gambro	Polyflux 210H	2.1	125	281	259	249	183	339	303	289	203	–	–	–	–	1500
Baxter	Exeltra Plus 210	2.1	125	287	277	252	202	350	328	292	222	–	–	–	–	1714

	Dialyzer	Surface Area (m2)	Priming Volume (ml)	Urea	Creatinine	Phosphate	Vit B12	Urea	Creatinine	Phosphate	Vit B12	Urea	Creatinine	Phosphate	Vit B12	KoA Urea
Fresenius	Optiflux F200NR	2.0	112	277	253	250	173	330	289	290	189	364	325	319	206	1317
Gambro	Polyflux 210H	2.1	125	281	259	249	183	339	303	289	203	–	–	–	–	1500
Baxter	Exeltra Plus 210	2.1	125	287	277	252	202	350	328	292	222	–	–	–	–	1714
Fresenius	Optiflux F250NR	2.5	135	289	273	273	206	–	–	–	–	383	357	356	251	1808