

Acid Concentrates

Acid concentrates can be produced for the following ranges:

Sodium 70-110 mEq/L

Potassium 0-5 mEq/L

Chloride 76.5-116mEq/L

Calcium 0-5 mEq/L

Magnesium 0-2 mEq/L

Acetic Acid 2-5 mEq/L

Dextrose 0-350 mg/dL



ACID CONCENTRATES								Corresponding Bicarbonate
Product Code	Na ⁺ mEq/L	Ca ⁺⁺ mEq/L	K ⁺ mEq/L	Mg ⁺⁺ mEq/L	Cl ⁻ mEq/L	Acetic Acid mEq/L	Dextrose mg/dL	
AC-100	81	2.5	2	1	86.5	4	100	BC-100
AC-101	81	2.5	2	1	86.5	4	200	BC-100
AC-102	80	3	2	1	86	4	200	BC-100
AC-103	80	2.5	2	1	85.5	4	200	BC-100
AC-104	81	2.5	0	1	84.5	4	200	BC-100
AC-105	80	3.5	2	1	86.5	4	200	BC-100
AC-106	79	3.5	2	1.5	86	4	200	BC-100
AC-107	80	2	2	1.5	85.5	4	200	BC-100
AC-109	80	2	3	1	86	4	200	BC-100
AC-110	80	3	3	1	87	4	200	BC-100
AC-111	80	3.5	0	1	84.5	4	200	BC-100
AC-112	80	3.5	1	1	85.5	4	200	BC-100
AC-113	80	2.5	4	1	87.5	4	200	BC-100
AC-114	80	2.5	3	1	86.5	4	200	BC-100
AC-115	81	2.5	1	1	85.5	4	200	BC-100
AC-116	81	3	2	1	87	4	200	BC-100
AC-200	100	2.5	2	1	105.5	4	100	BC-200
AC-201	100	2	2	1	105	4	100	BC-200
AC-202	100	2.5	1	1	104.5	4	100	BC-200
AC-203	100	2.5	3	1	106.5	4	100	BC-200
AC-204	100	3.5	2	1	106.5	4	100	BC-200
AC-205	100	2.5	2	0.75	105.25	4	200	BC-200
AC-206	100	3	2	0.75	105.75	4	200	BC-200
AC-207	100	2	2	0.75	104.75	4	200	BC-200
AC-208	100	3	2	0.75	105.75	2	100	BC-200
AC-209	100	3.5	2	0.75	106.25	2	100	BC-200
AC-210	100	3	1	0.75	105.75	2	100	BC-200
AC-211	100	3.5	2	0.75	106.25	4	200	BC-200
AC-212	100	2.5	3	0.75	106.25	4	100	BC-200
AC-213	100	0	0	0.75	100.75	4	200	BC-200
AC-214	100	2	1	1	104	4	200	BC-200
AC-215	100	2	3	0.75	105.75	4	200	BC-200
AC-216	100	2.5	3	0.75	106.25	4	200	BC-200
AC-217	100	3	1	0.75	104.75	4	200	BC-200
AC-218	100	3	3	0.75	106.75	4	200	BC-200
AC-219	100	2.25	0	1	103.25	4	100	BC-200
AC-220	100	0	1	1	102	4	100	BC-200
AC-221	100	2	1	1	104	4	100	BC-200
AC-222	100	2.25	1	1	104.25	4	100	BC-200
AC-223	100	3	1	1	105	4	100	BC-200
AC-224	100	2.25	2	1	105.25	4	100	BC-200
AC-225	100	3	2	1	106	4	100	BC-200
AC-226	100	2	3	1	106	4	100	BC-200
AC-227	100	2.25	3	1	106.25	4	100	BC-200
AC-228	100	3	3	1	107	4	100	BC-200
AC-229	100	2.25	4	1	107.25	4	100	BC-200
AC-300	103	2.5	1	1	107.5	3	200	BC-300
AC-301	105	2.5	2	1	110.5	3	200	BC-300
AC-302	103	2.5	2	1	108.5	3	100	BC-300
AC-303	103	3.5	3	1	110.5	3	100	BC-300
AC-304	103	3	2	1	109	3	100	BC-300
AC-305	103	3.5	2	1	109.5	3	100	BC-300
AC-306	103	3	1	1	108	3	100	BC-300
AC-307	103	2	2	0.75	107.75	3	200	BC-300
AC-308	105	3	2	1	111	4	150	BC-300
AC-309	103	2.5	2	0.75	108.25	3	200	BC-300
AC-310	103	2.5	3	1	109.5	3	200	BC-300
AC-311	103	3.5	3	1	110.5	3	200	BC-300