

Liquid Bicarbonate Concentrates



Liquid bicarbonate concentrates can be produced as:

Sodium 35.59 mEq/L

Bicarbonate 35.59 mEq/L

Chloride 0.20 mEq/L

BICARBONATE CONCENTRATES – LIQUID					Corresponding
Product Code	Dilution	Sodium	Chloride	Bicarbonate	Acids
BC-100	20.13X	59	20	39	36.83x
BC-201	26.14X	37	0	37	45x
BC-300	28.6X	35	0	35	35x