

## Blood Line Set



Medical grade PVC tubing minimizes side effects to health.

All components are manufactured and monitored consistently.

A variety of models are available for most brands of hemodialysis machines in operation worldwide; customized blood lines are available upon request.

Sterilization: E.O. gas and gamma ray.

CCC and CE certification.

BLOOD LINE SETS		
Code	Pump Segment (ID x OD x Length)mm	Descriptions
EDIM-BL-001	8.00 x 12.00 x 350	
EDIM-BL-002	6.3 x 9.8 x 400	
EDIM-BL-003	8.00 x 12.00 x 400	
EDIM-BL-004	8.00 x 12.00 x 400	
EDIM-BL-005	8.00 x 12.00 x 350	Big drip chamber
EDIM-BL-006	8.00 x 12.00 x 350	Small drip chamber
EDIM-BL-007	8.00 x 12.00 x 350	
EDIM-BL-008	8.00 x 12.00 x 350	