

Dry Bicarbonate



DRY BICARBONATE					
Product Code	Packets to Mix	Packaging	Specifications	Components Weight / Bag	Corresponding Acid
DB-100-2.5	2.5 Gallons	20 Bags/Case	Na 59mEq/L Cl 20mEq/L HCO ₃ 39m/Eq/L	NaCl 224g NaHCO ₃ 627g	36.83x
DB-100-15	15 Gallons	3 Bags/Case	Na 59mEq/L Cl 20mEq/L HCO ₃ 39m/Eq/L	NaCl 1336g NaHCO ₃ 3745g	36.83x
DB-100-25	25 Gallons	2 Bags/Case	Na 59mEq/L Cl 20mEq/L HCO ₃ 39m/Eq/L	NaCl 2227g NaHCO ₃ 6242g	36.83x
DB-201-2.1	2.1 Gallons	20 Bags/Case	Na 37mEq/L HCO ₃ 37m/Eq/L	NaHCO ₃ 650g	45x
DB-201-15	15 Gallons	3 Bags/Case	Na 37mEq/L HCO ₃ 37m/Eq/L	NaHCO ₃ 4617g	45x
DB-201-25	25 Gallons	2 Bags/Case	Na 37mEq/L HCO ₃ 37m/Eq/L	NaHCO ₃ 7696g	45x
DB-300-2.5	2.5 Gallons	20 Bags/Case	Na 35mEq/L HCO ₃ 35m/Eq/L	NaHCO ₃ 795g	35x
DB-300-15	15 Gallons	3 Bags/Case	Na 35mEq/L HCO ₃ 35m/Eq/L	NaHCO ₃ 4751g	35x
DB-300-25	25 Gallons	2 Bags/Case	Na 35mEq/L HCO ₃ 35m/Eq/L	NaHCO ₃ 7920g	35x