


A.V. Fistula Needle Set



Ultra thin wall and siliconized for smoother penetration.
 Special process on blade reduces pain and skin damage.
 Medical grade PVC tubing minimizes side effects to health.
 Sterilization: E.O.gas and gamma ray.
 CE certification.

DISPOSABLE A.V. FISTULA NEEDLE SETS					
Code	Spec.	Needle Length	Back Eye	Wing	Tubing
EDIM-AVF-001	16G	25±2.2mm	Y	Fixed	ID3.3xOD5.2x300mm
EDIM-AVF-002	16G	25±2.2mm	Y	Rotatable	ID3.3xOD5.2x300mm
EDIM-AVF-003	16G	32±2.2mm	Y	Fixed	ID3.3xOD5.2x300mm
EDIM-AVF-004	16G	32±2.2mm	Y	Rotatable	ID3.3xOD5.2x300mm
EDIM-AVF-005	17G	25±2.2mm	Y	Fixed	ID3.3xOD5.2x300mm
EDIM-AVF-006	17G	25±2.2mm	Y	Rotatable	ID3.3xOD5.2x300mm
EDIM-AVF-007	17G	32±2.2mm	Y	Fixed	ID3.3xOD5.2x300mm
EDIM-AVF-008	17G	32±2.2mm	Y	Rotatable	ID3.3xOD5.2x300mm
EDIM-AVF-009	15G	25±2.2mm	Y	Fixed	ID3.3xOD5.2x300mm
EDIM-AVF-010	15G	25±2.2mm	Y	Rotatable	ID3.3xOD5.2x300mm
EDIM-AVF-011	15G	32±2.2mm	Y	Fixed	ID3.3xOD5.2x300mm
EDIM-AVF-012	15G	32±2.2mm	Y	Rotatable	ID3.3xOD5.2x300mm